

ARENA™

FIXATION BUTTON CONTINUOUS LOOP





Product Information

ARENA - FIXATION BUTTON CONTINUOUS LOOP

ARENA is a cortical fixation button ideal for either anatomical or traditional Anterior Cruciate Ligament (ACL) reconstruction. With its titanium button and ultra high molecular weight polyethylene (UHMWPE) continuous loop, ARENA offers a safe and strong ligament fixation with an ultimate load failure 50% higher than the market leader.*

ARENA is preloaded with one white and one coloured suture to both lead and flip the button, sizes range from 12mm to 60mm.




*Data on file

PRODUCT CODES

ARENA™

Product Code	Description
D201002	ARENA 12mm
D201003	ARENA 15mm
D201004	ARENA 20mm
D201005	ARENA 25mm
D201006	ARENA 30mm
D201007	ARENA 35mm
D201008	ARENA 40mm
D201009	ARENA 45mm
D201010	ARENA 50mm
D201011	ARENA 55mm
D201012	ARENA 60mm



DISTRIBUTED BY



Head Office
 Industrial Park
 Bourton-on-the-Water
 Gloucestershire
 GL54 2HQ
 United Kingdom

tel. +44 (0) 1451 821311
 email info@orthod.com
 website www.orthod.com