

# HIVAC™ 7

CEMENTATION



## Product Information

## HIVAC7 - BONE CEMENT MIXING AND DELIVERY SYSTEM

HIVAC7 is a bone cement mixing and delivery device that uses a vacuum level of 550mmHg to ensure a good mix and optimal bone cement porosity<sup>1</sup>.

### BENEFITS OF HIVAC7

#### SIMPLE AND EASY TO USE

HIVAC7 has few components and is simple and easy to use, allowing theatre staff to feel confident in producing a high quality mix every time.

#### INTUITIVE DELIVERY GUN

HIVAC7 can easily be inserted into the delivery gun. The rod only has teeth on one side and can easily be engaged or disengaged.

#### SNAP-OFF POINT

The mixing rod slides smoothly within the mixing chamber and features a coloured indicator to allow the snap point to be clearly visible.

#### CLEAR MIXING TUBE

HIVAC7 features a clear tube permitting an unclouded view of the cement during mixing. This offers peace of mind to the user as mix quality can be visually assessed prior to delivery.

#### TRUSTED PRODUCTS

Our cementation range has been used in millions of procedures for more than 20 years across more than 40 countries.

#### FEMORAL PRESSURISER

Flexible material achieves a seal on a wide range of femurs to allow pressurisation. Easy to fit onto HIVAC7.



### PRODUCT CODES

#### HIVAC™7

Product Code	Description	Box Qty
C701	HIVAC7 (80g capacity; 40g Simplex)	10
CP701	HIVAC7 with Pressuriser	10
C702	HIVAC7 (Qty 2 in pack)	10
C703	HIVAC7 Triple Mix (120g capacity; 80g Simplex)	5
C800	HIVAC7 Delivery Syringe	20
CB900	HIVAC7 Combination Hip Pack	5
H550	High Vacuum Pump	1
H759	Ratchet Delivery Gun for HIVAC7	1
H759L	Ratchet Delivery Gun for HIVAC7 Triple Mix	1
P721	Femoral Cement Pressuriser	10

REF 1 Dunne N.J, Orr J.F, Mushipe M.T, Eveleigh R.J, 2002. The relationship between porosity and fatigue characteristics of bone cement. *Biomaterials* 24 (2003) 238-245.

#### DISTRIBUTED BY



Head Office  
 Industrial Park  
 Bourton-on-the-Water  
 Gloucestershire  
 GL54 2HQ  
 United Kingdom

tel. +44 (0) 1451 821311  
 email [info@orthod.com](mailto:info@orthod.com)  
 website [www.orthod.com](http://www.orthod.com)